

C-C motif chemokine 5

Catalog No: #AP81496



Package Size: #AP81496-1 50ug #AP81496-2 100ug #AP81496-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	SIS-delta, Small-inducible cytokine A5, T-cell-specific protein RANTES, SCYA5
Accession No.	Uniprot ID:O97919
Uniprot	O97919
GeneID	327712;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CCL5

Note: This product is for in vitro research use only