

Hematopoietic progenitor cell antigen CD34

Catalog No: #AP81498



Package Size: #AP81498-1 50ug #AP81498-2 100ug #AP81498-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hematopoietic progenitor cell antigen CD34
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CD34
Accession No.	Uniprot ID:Q64314
Uniprot	Q64314
GeneID	12490;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cd34

Note: This product is for in vitro research use only