

Adenosine deaminase CECR1

Catalog No: #AP81506



Package Size: #AP81506-1 50ug #AP81506-2 100ug #AP81506-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenosine deaminase CECR1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cat eye syndrome critical region protein 1,ADA2,ADGF,IDGFL
Accession No.	Uniprot ID:Q9NZK5
Uniprot	Q9NZK5
GeneID	51816;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CECR1

Note: This product is for in vitro research use only