

# G1/S-specific cyclin-D1

Catalog No: #AP81508



Package Size: #AP81508-1 50ug #AP81508-2 100ug #AP81508-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	G1/S-specific cyclin-D1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Accession No.	Uniprot ID:Q64HP0
Uniprot	Q64HP0
GeneID	449028;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CCND1

Note: This product is for in vitro research use only