

C-C chemokine receptor type 2

Catalog No: #AP81509



Package Size: #AP81509-1 50ug #AP81509-2 100ug #AP81509-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C chemokine receptor type 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C-C CKR-2, Monocyte chemoattractant protein 1 receptor, MCP-1-R, CD192, CMKBR2
Accession No.	Uniprot ID:P41597
Uniprot	P41597
GeneID	729230;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CCR2

Note: This product is for in vitro research use only