

C-C chemokine receptor type 6

Catalog No: #AP81514



Package Size: #AP81514-1 50ug #AP81514-2 100ug #AP81514-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C chemokine receptor type 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C-C CKR-6,KY411,CD196,Cmkbr6
Accession No.	Uniprot ID:O54689
Uniprot	O54689
GeneID	12458;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccr6

Note: This product is for in vitro research use only