

Cyclin-dependent kinase 4

Catalog No: #AP81518



Package Size: #AP81518-1 50ug #AP81518-2 100ug #AP81518-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Cell division protein kinase 4
Accession No.	Uniprot ID:P79432
Uniprot	P79432
GeneID	100144492;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CDK4

Note: This product is for in vitro research use only