

Carbonic anhydrase 2

Catalog No: #AP81519



Package Size: #AP81519-1 50ug #AP81519-2 100ug #AP81519-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carbonic anhydrase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Carbonate dehydratase II,Carbonic anhydrase II,CA-II,Car2
Accession No.	Uniprot ID:P00920
Uniprot	P00920
GeneID	12349;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ca2

Note: This product is for in vitro research use only