

Calreticulin

Catalog No: #AP81522



Package Size: #AP81522-1 50ug #AP81522-2 100ug #AP81522-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calreticulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CRP55, Calregulin, Endoplasmic reticulum resident protein 60, ERp60, HACBP
Accession No.	Uniprot ID:P28491
Uniprot	P28491
GeneID	100381266;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CALR

Note: This product is for in vitro research use only