

Macrophage-capping protein

Catalog No: #AP81523



Package Size: #AP81523-1 50ug #AP81523-2 100ug #AP81523-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Macrophage-capping protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Actin regulatory protein CAP-G,Actin-capping protein GCAP39,Myc basic motif homolog 1,Mbh1
Accession No.	Uniprot ID:P24452
Uniprot	P24452
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Capg

Note: This product is for in vitro research use only