

# CD276 antigen

Catalog No: #AP81524



Package Size: #AP81524-1 50ug #AP81524-2 100ug #AP81524-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD276 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	4lg-B7-H3,B7 homolog 3,B7-H3,Costimulatory molecule,CD276,B7H3,PSEC0249,UNQ309/PRO352
Accession No.	Uniprot ID:Q5ZPR3
Uniprot	Q5ZPR3
GeneID	80381;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD276

Note: This product is for in vitro research use only