

Cyclin-dependent kinase 5 activator 1

Catalog No: #AP81525



Package Size: #AP81525-1 50ug #AP81525-2 100ug #AP81525-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 5 activator 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CDK5 activator 1,Cyclin-dependent kinase 5 regulatory subunit 1,TPKII regulatory subunit,Cdk5r,Nck5a
Accession No.	Uniprot ID:P61809
Uniprot	P61809
GeneID	12569;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdk5r1

Note: This product is for in vitro research use only