

## C-C motif chemokine 3

Catalog No: #AP81527

Package Size: #AP81527-1 50ug #AP81527-2 100ug #AP81527-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-C motif chemokine 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Macrophage inflammatory protein 1-alpha,MIP-1-alpha,Small-inducible cytokine A3,SCYA3
Accession No.	Uniprot ID:Q8SQA6
Uniprot	Q8SQA6
GeneID	282170;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CCL3

Note: This product is for in vitro research use only