

Platelet glycoprotein 4

Catalog No: #AP81530



Package Size: #AP81530-1 50ug #AP81530-2 100ug #AP81530-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet glycoprotein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fatty acid translocase,FAT,Glycoprotein IIIb,GPIIIB,Leukocyte differentiation antigen CD36,PAS IV,PAS-4,Platelet collagen receptor,Platelet glycoprotein IV,GPIV,Thrombospondin receptor,CD36,GP3B,GP4
Accession No.	Uniprot ID:P16671
Uniprot	P16671
GeneID	948;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD36

Note: This product is for in vitro research use only