

## CD5 antigen-like

Catalog No: #AP81533

Package Size: #AP81533-1 50ug #AP81533-2 100ug #AP81533-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD5 antigen-like
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Apoptosis inhibitor expressed by macrophages,hAIM,CT-2,IgM-associated peptide,SP-alpha,API6,UNQ203/PRO229
Accession No.	Uniprot ID:O43866
Uniprot	O43866
GeneID	922;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD5L

Note: This product is for in vitro research use only