

Cyclin-dependent-like kinase 5

Catalog No: #AP81535



Package Size: #AP81535-1 50ug #AP81535-2 100ug #AP81535-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent-like kinase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CR6 protein kinase,CRK6,Cell division protein kinase 5,Serine/threonine-protein kinase PSSALRE,Tau protein kinase II catalytic subunit,TPKII catalytic subunit,Cdkn5,Crk6
Accession No.	Uniprot ID:P49615
Uniprot	P49615
GeneID	12568;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdk5

Note: This product is for in vitro research use only