

CD40 ligand

Catalog No: #AP81536



Package Size: #AP81536-1 50ug #AP81536-2 100ug #AP81536-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD40 ligand
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CD40-L, Tumor necrosis factor ligand superfamily member 5, CD154, CD40L, TNFSF5
Accession No.	Uniprot ID:P51749
Uniprot	P51749
GeneID	282387;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CD40LG

Note: This product is for in vitro research use only