

T-cell surface glycoprotein CD8 beta chain

Catalog No: #AP81538



Package Size: #AP81538-1 50ug #AP81538-2 100ug #AP81538-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD8 beta chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CD8 antigen 37 kDa chain,OX-8 membrane antigen,CD8b,Cd8b1
Accession No.	Uniprot ID:P05541
Uniprot	P05541
GeneID	24931;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd8b

Note: This product is for in vitro research use only