

Cyclin-dependent kinase 1

Catalog No: #AP81539

Package Size: #AP81539-1 50ug #AP81539-2 100ug #AP81539-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CDK1,Cell division control protein 2 homolog,Cell division protein kinase 1,p34 protein kinase,CDC2,CDC28A,CDKN1,P34CDC2
Accession No.	Uniprot ID:P06493
Uniprot	P06493
GeneID	983;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDK1

Note: This product is for in vitro research use only