

# CD40 ligand

Catalog No: #AP81540



Package Size: #AP81540-1 50ug #AP81540-2 100ug #AP81540-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD40 ligand
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	CD40-L, Tumor necrosis factor ligand superfamily member 5, CD40L, TNFSF5, CD154
Accession No.	Uniprot ID:O97626
Uniprot	O97626
GeneID	403468;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: CD40LG

Note: This product is for in vitro research use only