

## Cyclin-dependent kinase 4 inhibitor C

Catalog No: #AP81541



Package Size: #AP81541-1 50ug #AP81541-2 100ug #AP81541-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cyclin-dependent kinase 4 inhibitor C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cyclin-dependent kinase 6 inhibitor,p18-INK4c,p18-INK6,CDKN6
Accession No.	Uniprot ID:P42773
Uniprot	P42773
GeneID	1031;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CDKN2C

Note: This product is for in vitro research use only