

Cyclin-dependent kinase 7

Catalog No: #AP81544

Package Size: #AP81544-1 50ug #AP81544-2 100ug #AP81544-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	39 kDa protein kinase,p39 Mo15,CDK-activating kinase 1,Cell division protein kinase 7,Serine/threonine-protein kinase 1,TFIIH basal transcription factor complex kinase subunit,CAK,CAK1,CDKN7,MO15,STK1
Accession No.	Uniprot ID:P50613
Uniprot	P50613
GeneID	1022;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDK7

Note: This product is for in vitro research use only