

Cell division control protein 42 homolog

Catalog No: #AP81546



Package Size: #AP81546-1 50ug #AP81546-2 100ug #AP81546-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cell division control protein 42 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	G25K GTP-binding protein
Accession No.	Uniprot ID:Q90694
Uniprot	Q90694
GeneID	395917;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CDC42

Note: This product is for in vitro research use only