

Baculoviral IAP repeat-containing protein 1a

Catalog No: #AP81547



Package Size: #AP81547-1 50ug #AP81547-2 100ug #AP81547-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Baculoviral IAP repeat-containing protein 1a |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | Neuronal apoptosis inhibitory protein 1,Birc1a,Naip |
| Accession No. | Uniprot ID:Q9QWK5 |
| Uniprot | Q9QWK5 |
| GeneID | 17940; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Naip1

Note: This product is for in vitro research use only