

Complement factor H

Catalog No: #AP81550

Package Size: #AP81550-1 50ug #AP81550-2 100ug #AP81550-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement factor H
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein beta-1-H,Hf1
Accession No.	Uniprot ID:P06909
Uniprot	P06909
GeneID	12628;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cfh

Note: This product is for in vitro research use only