

Major centromere autoantigen B

Catalog No: #AP81551

Package Size: #AP81551-1 50ug #AP81551-2 100ug #AP81551-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Major centromere autoantigen B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Centromere protein B,CENP-B,Cenp-b
Accession No.	Uniprot ID:P27790
Uniprot	P27790
GeneID	12616;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cenpb

Note: This product is for in vitro research use only