

C-C chemokine receptor type 7

Catalog No: #AP81553



Package Size: #AP81553-1 50ug #AP81553-2 100ug #AP81553-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C chemokine receptor type 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C-C CKR-7, Epstein-Barr virus-induced G-protein coupled receptor 1, EB11, MIP-3 beta receptor, CD197, Cmkbr7, Ebi1, Ebi1h
Accession No.	Uniprot ID:P47774
Uniprot	P47774
GeneID	12775;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ccr7

Note: This product is for in vitro research use only