

# T-cell surface antigen CD2

Catalog No: #AP81554



Package Size: #AP81554-1 50ug #AP81554-2 100ug #AP81554-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	T-cell surface antigen CD2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LFA-2,LFA-3 receptor,OX-34 antigen,T-cell surface antigen T11/Leu-5,CD2
Accession No.	Uniprot ID:P08921
Uniprot	P08921
GeneID	497761;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cd2

Note: This product is for in vitro research use only