

T-cell surface glycoprotein CD3 epsilon chain

Catalog No: #AP81555



Package Size: #AP81555-1 50ug #AP81555-2 100ug #AP81555-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-cell surface glycoprotein CD3 epsilon chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CD3e
Accession No.	Uniprot ID:Q28073
Uniprot	Q28073
GeneID	281054;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CD3E

Note: This product is for in vitro research use only