

Calpain-2 catalytic subunit

Catalog No: #AP81559



Package Size: #AP81559-1 50ug #AP81559-2 100ug #AP81559-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calpain-2 catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	80 kDa M-calpain subunit,CALP80,Calcium-activated neutral proteinase 2,CANP 2,Calpain M-type,Calpain-2 large subunit,Millimolar-calpain,M-calpain
Accession No.	Uniprot ID:O08529
Uniprot	O08529
GeneID	12334;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Capn2

Note: This product is for in vitro research use only