

Calpain-5

Catalog No: #AP81560



Package Size: #AP81560-1 50ug #AP81560-2 100ug #AP81560-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Calpain-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	New calpain 3,nCL-3,Ncl3
Accession No.	Uniprot ID:O08688
Uniprot	O08688
GeneID	12337;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Capn5

Note: This product is for in vitro research use only