

Early activation antigen CD69

Catalog No: #AP81561



Package Size: #AP81561-1 50ug #AP81561-2 100ug #AP81561-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Early activation antigen CD69
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Activation inducer molecule,AIM,BL-AC/P26,C-type lectin domain family 2 member C,EA1,Early T-cell activation antigen p60,GP32/28,Leukocyte surface antigen Leu-23,MLR-3,CD69,CLEC2C
Accession No.	Uniprot ID:Q07108
Uniprot	Q07108
GeneID	969;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD69

Note: This product is for in vitro research use only