

## T-cell surface glycoprotein CD8 beta chain

Catalog No: #AP81564



Package Size: #AP81564-1 50ug #AP81564-2 100ug #AP81564-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	T-cell surface glycoprotein CD8 beta chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lymphocyte antigen 3,T-cell membrane glycoprotein Ly-3,T-cell surface glycoprotein Lyt-3,CD8b,Cd8b1,Ly-3,Lyt-3,Lyt3
Accession No.	Uniprot ID:P10300
Uniprot	P10300
GeneID	12526;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cd8b

Note: This product is for in vitro research use only