

Protein chibby homolog 1

Catalog No: #AP81566



Package Size: #AP81566-1 50ug #AP81566-2 100ug #AP81566-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein chibby homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytosolic leucine-rich protein,PIGEA-14,PKD2 interactor, Golgi and endoplasmic reticulum-associated 1,Cby,Pgea1
Accession No.	Uniprot ID:Q8K4I6
Uniprot	Q8K4I6
GeneID	246768;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cby1

Note: This product is for in vitro research use only