

C-C motif chemokine 21a

Catalog No: #AP81567



Package Size: #AP81567-1 50ug #AP81567-2 100ug #AP81567-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 21a
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	6Ckine,Beta-chemokine exodus-2,Small-inducible cytokine A21a,Thymus-derived chemotactic agent 4,TCA4,Scya21,Scya21a
Accession No.	Uniprot ID:P84444
Uniprot	P84444
GeneID	100504362;18829;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccl21a

Note: This product is for in vitro research use only