

C-C motif chemokine 18

Catalog No: #AP81568



Package Size: #AP81568-1 50ug #AP81568-2 100ug #AP81568-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Alternative macrophage activation-associated CC chemokine 1,AMAC-1,CC chemokine PARC,Dendritic cell chemokine 1,DC-CK1,Macrophage inflammatory protein 4,MIP-4,Pulmonary and activation-regulated chemokine,Small-inducible cytokine A18,AMAC1,DCCK1,MIP4,PARC,SCYA18
Accession No.	Uniprot ID:P55774
Uniprot	P55774
GeneID	6362;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCL18

Note: This product is for in vitro research use only