

CD166 antigen

Catalog No: #AP81570



Package Size: #AP81570-1 50ug #AP81570-2 100ug #AP81570-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD166 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Activated leukocyte cell adhesion molecule,BEN glycoprotein,Protein DM-GRASP,Protein JC7,SC1 glycoprotein,CD166
Accession No.	Uniprot ID:P42292/ncbi
Uniprot	P42292
GeneID	396092;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ALCAM

Note: This product is for in vitro research use only