

Platelet glycoprotein 4

Catalog No: #AP81571



Package Size: #AP81571-1 50ug #AP81571-2 100ug #AP81571-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet glycoprotein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Adipocyte membrane protein,Fatty acid translocase,Fatty acid transport protein,Glycoprotein IIIb,GPIIIB,PAS IV,PAS-4,Platelet glycoprotein IV,GPIV,CD36,Fat
Accession No.	Uniprot ID:Q07969
Uniprot	Q07969
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd36

Note: This product is for in vitro research use only