

CD40 ligand

Catalog No: #AP81572



Package Size: #AP81572-1 50ug #AP81572-2 100ug #AP81572-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD40 ligand
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CD40-L, Tumor necrosis factor ligand superfamily member 5, CD154, Cd40l, Tnfsf5
Accession No.	Uniprot ID:Q9Z2V2
Uniprot	Q9Z2V2
GeneID	84349;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cd40lg

Note: This product is for in vitro research use only