

Chromogranin-A

Catalog No: #AP81574

Package Size: #AP81574-1 50ug #AP81574-2 100ug #AP81574-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chromogranin-A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CgA,Pituitary secretory protein I,SP-I
Accession No.	Uniprot ID:P05059
Uniprot	P05059
GeneID	281070;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CHGA

Note: This product is for in vitro research use only