

C-C motif chemokine 23

Catalog No: #AP81577

Package Size: #AP81577-1 50ug #AP81577-2 100ug #AP81577-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 23
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CK-beta-8,CKB-8,Macrophage inflammatory protein 3,MIP-3,Myeloid progenitor inhibitory factor 1,MPIF-1,Small-inducible cytokine A23,MIP3,MPIF1,SCYA23
Accession No.	Uniprot ID:P55773
Uniprot	P55773
GeneID	6368;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CCL23

Note: This product is for in vitro research use only