

## CAMPATH-1 antigen

Catalog No: #AP81579



Package Size: #AP81579-1 50ug #AP81579-2 100ug #AP81579-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CAMPATH-1 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CDw52,Cambridge pathology 1 antigen,Epididymal secretory protein E5,Human epididymis-specific protein 5,He5,CD52,CDW52,HE5
Accession No.	Uniprot ID:P31358
Uniprot	P31358
GeneID	1043;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD52

Note: This product is for in vitro research use only