

Calpain-3

Catalog No: #AP81584



Package Size: #AP81584-1 50ug #AP81584-2 100ug #AP81584-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calpain-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Calcium-activated neutral proteinase 3,CANP 3,Calpain L3,Calpain p94,Muscle-specific calcium-activated neutral protease 3,New calpain 1,nCL-1,CANP3,CANPL3,NCL1
Accession No.	Uniprot ID:P20807
Uniprot	P20807
GeneID	825;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CAPN3

Note: This product is for in vitro research use only