

C-C motif chemokine 25

Catalog No: #AP81586



Package Size: #AP81586-1 50ug #AP81586-2 100ug #AP81586-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 25
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Chemokine TECK, Small-inducible cytokine A25, Thymus-expressed chemokine, SCYA25, TECK
Accession No.	Uniprot ID:O15444
Uniprot	O15444
GeneID	6370;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CCL25

Note: This product is for in vitro research use only