

Cyclin-dependent kinase 7

Catalog No: #AP81590



Package Size: #AP81590-1 50ug #AP81590-2 100ug #AP81590-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	39 kDa protein kinase,P39 Mo15,CDK-activating kinase,CR4 protein kinase,CRK4,Cell division protein kinase 7,Protein-tyrosine kinase MPK-7,TFIIH basal transcription factor complex kinase subunit,Cak,Cdkn7,Crk4,Mo15,Mpk-7
Accession No.	Uniprot ID:Q03147
Uniprot	Q03147
GeneID	12572;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdk7

Note: This product is for in vitro research use only