

## Cysteine dioxygenase type 1

Catalog No: #AP81591



Package Size: #AP81591-1 50ug #AP81591-2 100ug #AP81591-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cysteine dioxygenase type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Cysteine dioxygenase type I,CDO
Accession No.	Uniprot ID:Q3SZU4
Uniprot	Q3SZU4
GeneID	514462;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CDO1

Note: This product is for in vitro research use only