

Homeobox protein CDX-2

Catalog No: #AP81592



Package Size: #AP81592-1 50ug #AP81592-2 100ug #AP81592-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Homeobox protein CDX-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CDX-3,Caudal-type homeobox protein 2,CDX3
Accession No.	Uniprot ID:Q99626
Uniprot	Q99626
GeneID	1045;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDX2

Note: This product is for in vitro research use only