

CD5 antigen-like

Catalog No: #AP81594



Package Size: #AP81594-1 50ug #AP81594-2 100ug #AP81594-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD5 antigen-like
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Apoptosis inhibitor expressed by macrophages,mAIM,Apoptosis inhibitory 6,SP-alpha,Aim,Api6
Accession No.	Uniprot ID:Q9QWK4
Uniprot	Q9QWK4
GeneID	11801;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd5l

Note: This product is for in vitro research use only