

C-C motif chemokine 27

Catalog No: #AP81597

Package Size: #AP81597-1 50ug #AP81597-2 100ug #AP81597-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 27
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CC chemokine ILC,Cutaneous T-cell-attracting chemokine,CTACK,ESkine,IL-11 R-alpha-locus chemokine,ALP,Skinkine,Small-inducible cytokine A27,Ilc,Scya27
Accession No.	Uniprot ID:Q9Z1X0
Uniprot	Q9Z1X0
GeneID	100039863;100039939;100040048;100861978;20301;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccl27

Note: This product is for in vitro research use only