

Cyclin-dependent kinase 11A

Catalog No: #AP81598



Package Size: #AP81598-1 50ug #AP81598-2 100ug #AP81598-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Cyclin-dependent kinase 11A |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Cell division cycle 2-like protein kinase 2, Cell division protein kinase 11A, Galactosyltransferase-associated protein kinase p58/GTA, PITSLRE serine/threonine-protein kinase CDC2L2, CDC2L2, CDC2L3, PITSLREB |
| Accession No. | Uniprot ID:Q9UQ88 |
| Uniprot | Q9UQ88 |
| GeneID | 728642; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CDK11A

Note: This product is for in vitro research use only